Instruction manual
Manual de Instrucciones

Master-care®

The Swedish Back Care System™

Back-A-Traction®

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(CN-A1, CN-A2, CN-A3, CN-AH6 and/y CN-B1, CN-B2, CN-B3)

Manufactured under license for Swedish Backcare System, Inc.
Bajo licencia de manufactura para Swedish Backcare System, Inc.

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ENGLISH

We congratulate you on your decision to seek a healthier life through the use of Back-A-Traction from Mastercare.

Mastercare Back-A-Traction is a gravity traction table that will help you relieve and prevent back pain and the wearing out of joints. It is a good tool to help you regain your back health and to maintain it for the rest of your life.

Mastercare Back-A-Traction's unique feature, the Sliding Backrest, moves on rollers on the steel frame. This patented feature makes it possible to obtain optimal beneficial traction at a low inversion angle. While your spine is extended, the back is perfectly supported allowing the muscles to relax completely. You will experience an unloading of the pressure of gravitational forces from your discs and joints.

To get optimal results, you should alternate stretching exercises and total relaxation. The passive traction and the stretching exercises on the table involve stretching of the muscles, muscular attachments, tendons and crucial ligaments of your entire back and legs. In the beginning this may cause a slight stiffness, which is a sign that the mentioned parts of your body are not used to this kind of stretching.

Even if you become sore in the beginning, you may continue your exercises with less intensity, as the soreness should disappear as the muscles and tendons increase their flexibility. Daily use of Mastercare Back-A-Traction will lead to a healthier life.

CAUTION ADVISED IN CERTAIN MEDICAL PRE-EXISTING CONDITIONS

ALWAYS consult your doctor before starting a regular use of the bench if you are uncertain about the state of your health, or suffer from any of the following conditions:

Glaucoma, detached retina, conjunctivitis, uncontrolled hypertension (high blood pressure), impaired pulmonary functions, arrhythmia, carotid artery stenosis, severe vascular disease, aneurysm, tendency towards fainting or dizzy spells, cerebral sclerosis, recent stroke or transient ischemic attack, history of brain lesions, other heart or circulatory disorders for which you are being treated, severe back problems or spinal injuries, bone weakness (osteoporosis), recent unhealed fractures, medullary pins, artificial hip joints or other orthopedic supports surgically implanted, middle ear infection, recently had surgery, using special medication, ventral or hiatal hernia, tendency towards bleeding, extreme obesity, invertaphobia (fear of inversion itself), and pregnancy.

IMPORTANT SAFEGUARDS

Not for use by minors under the age of 18.

Have others remain clear while the table is in use, to avoid bodily injury from moving parts.

As with any therapy or exercise equipment, you should consult your doctor before regular use of Mastercare Back-A-Traction.

The table is designed for use by one person at a time, with both feet fastened in the foot-locking device. In the beginning, use Mastercare Back-A-Traction in the presence of a person who is familiar with its functions.

Before using Mastercare Back-A-Traction on your own, you must be certain that you fully understand how the table works.

Be sensitive to the signals from your body. Consult your doctor or physical therapist if you feel uncomfortable, have pain or notice any unusual changes that you think might prevent you from using Mastercare Back-A-Traction.

Do not use Mastercare Back-A-Traction immediately after a meal, or if you have a temperature, a cold, or severe head- or ear ache.

Mastercare Back-A-Traction is designed for persons with a height between 4,6 and 7 feet (135 to 215 cm), and a maximum weight of 330 lbs. (150 kg) (Exception; Model CN-AH6 Heavy Duty for Athletes and Sports Teams where maximum weight is 440lbs./200 kg)

Inspect your Mastercare unit and ALWAYS make sure that:

- the wires on the stand are stretched
- it stands on a level, horizontal surface
- all nuts and bolts are securely fastened
- the “fork” of the foot locking device is working properly/perfectly

Learn how the two thumb activated safety catches work and make sure that they are both engaged. If you temporarily disengage one of the safety catches when adjusting the height, you must re-engage the safety catch when you are done with the height adjustment. If both safety catches are disengaged, the backrest will swing without safety stops and you are “out of control” and may get injured.

When using Mastercare Back-A-Traction in prone position/laying face down (**ONLY pertaining Models A, NOT B**); you should always be assisted by another person, who is familiar with its operation.

Before stepping up on the foot plate, always make sure that the catch for the height adjustment is in a locked position.

If you do not balance the height adjustment well when operating Mastercare Back-A-Traction, re-adjust the height on the foot rest to obtain the best possible balance. (In order to make sure that the backrest will not rock back and forth when you do the stretching exercises, it is usually advantageous to have it adjusted so that the part where your upper body is resting weights a little more). Wiggle the backrest slightly when you push down the “push-buttons on the handles/hand grips to loosen the thumb activated safety catches. If the “push-buttons” works stiffly, the height is not properly adjusted. To correct this, the foot rest has to be extended or shortened.

If your back is in a very bad condition when you start using Mastercare Back-A-Traction, it will for the first couple of weeks, probably be enough for you to just stretch and relax in 15 degree with your arms over your head. Then, as your condition improves, you can begin to do the stretching exercise program.

A BASIC STRETCHING EXERCISE PROGRAM

Follow the instructions you have received and **watch the enclosed DVD for further understanding**. Be cautious and sensitive to the signals from your body. For best result you should relax and take some deep breaths between the different stretching exercises. It is also important that you relax as much as possible during each stretching exercise.

The following instruction will tell you to repeat certain exercises. It is for you to decide if you want to repeat them, and how many times. An easy start is recommended. Increase muscular stretching and relaxation with each repetition. Most of the following stretching exercises should be done, if possible, with your arms over your head.

1. You must first adjust the height for correct balance of the back rest according to your height. Press down one of the “push-buttons on the handles/hand grips (preferably the one on the opposite side) to release the safety catch and lock it in an “open” position by swinging the safety catch down on the side of the stand. This allows all adjustments to be made from one side.

2. Next, press down the other push-button to release the lock and incline/rotate the back rest to the horizontal position where it locks in place once you have released the “push-button.
3. To adjust the balance, pull the **height and balance adjustment knob** and slide the square-tube/bar according to the guideline markings provided. Test the back rest and fine tune the balance to your height and weight as necessary.

4. Next, press the “push-button” and rotate the back rest towards its upright position. It will lock at a slight incline. Now, return the locked safety catch to the operating position by pressing the push-button and releasing the safety catch.

You're now ready to use Mastercare Back-A-Traction. As a first time user you should always have someone help you until you feel comfortable using Mastercare Back-A-Traction.

5. Remove your shoes, and then unlock the foot restraint/foot-locking device, swinging the handle to the wide open position. Step into Mastercare Back-A-Traction, sitting against the back rest while placing your feet between the

restraining pads (foot-rolls). Secure your feet with the foot restraint, moving along its ratchet until your feet are held, comfortably and securely between the pads (foot-rolls).

6. Now, lean back onto the back rest and take hold of the handles/hand grips and press down the “push-buttons” with your thumbs. Thus releasing the thumb locks and then push back to the horizontal position where the back rest will lock in place once you have released the “push-button.”

Before exercising, you should stretch by pushing against the handles/hand grips. Throughout the session, remember to relax and breath deeply and allow Mastercare Back-A-Traction to work, stretching the tension from your muscles, tendons and ligaments, the common causes of back pain.

7. After stretching, press and release the “push-buttons” and incline/rotate further backwards to the first catch at 15 degrees, the position you should use for your first few weeks with Mastercare Back-A-Traction. (15 degree is sufficient and the degree we recommend to stay in – never need to go further – thanks to its moving back rest). Repeat the stretching procedure by pushing, as before, against the handles. Now, raise your hands slowly over your head and prepare to exercise. First begin with a few repetitions doing dorsiflexion with your feet (bending your feet towards you - so the toes are brought closer to the shin) under/in the restraint pads (foot-rolls).
8. Next, alternate bending your knees in a slow and easy walking motion. Each exercise should be repeated a few times for best effect. If difficult; bend both knees at the same time and then down with one knee at a time.

Remember to relax and breathe deeply between each exercise.

9. The exercise involves arching and lifting your buttocks off of the back rest. (Tighten your buttock, if possible; lift (arch) the buttock off the back rest – up and down).
10. Next, push your neck flat to the back rest while lifting/arching your shoulders upward.
11. After a few repetitions, place one hand on your stomach and then reach over your head with the other, as if trying to pick apples just out of reach.

Now, repeat the process, this time holding your arm out to your side. You should alternate your palms up and down from time to time during this exercise. Finally, repeat the exercise holding your hand near your waist.

From this finishing position, move right into the next exercise, moving your arm in circles as if swimming the backstroke. Repeat these hand and arm exercises with your other arm. End with your back flat to the backrest and both hands together on your stomach. (You can do both arms at the same time, if you wish)

12. The next exercise involves a three part move. Raise both hands over your head and then twist your torso and hips to the right, raising the right knee at the same time. Return to the flat back, flat knee position and twist to the left side. Return to the flat position and repeat the move a few times. End the exercise by bringing your hands back to rest on your stomach.
13. The final exercise involves quarter sit-ups. Lift your head with the chin towards your chest, leaving the rest of your back flat to the backrest. Proceed only as you feel comfortable to do so. If difficult; start only by nodding or pulling your chin in towards your neck down without lifting the head.

End by returning flat to the backrest and raising your arms over your head for a last reaching stretch.

14. Begin the dismount by taking hold of the handles/hand grips, push down the “push-buttons” to release the safety catch locks, and returning to the horizontal position. Relax for at least 60 seconds and move your knees in a walking mode for a few repetitions as before. Once again; push down the “push-buttons” to release the safety catch locks and return to the 15 degree upright position. Pause again for a moment.

15. Finally, push down the “push-buttons” to release the safety catch locks and return the back rest to its last and most upright locking point. You should be able to slip a foot out from under the foot pads (foot-rolls) and then release the foot restraint/foot locking device to its open position.

Step out and away from your Mastercare Back-A-Traction.

It is very important to relax after each session. Avoid lifting and strenuous movements, and take it easy for a few minutes.

Now you have a basic stretching exercise program.

Start with caution. In the beginning "hold" means a few seconds. After you have been training for some weeks it may mean 20 to 30 seconds. Notice the reactions of your body. Intensify the exercises gradually.

Feel free to talk to your back specialist about an exercise program just for you.

After some time, when you are familiar with Mastercare Back-A-Traction and the reactions of your body, you can do the stretching exercises in the 30 degree inversion angle. (The second position after horizontal). However, 15 degree has proven to be most effective and therefore all you need.

Total relaxation is necessary to get the best results.

ADDITIONAL STRETCHING EXERCISES

- A. Put your hands on your shoulders. Rotate arms and shoulders slowly forwards and backwards. Relax. Breathe. Repeat.
- B. Clasp one wrist and pull the arm over your head in the long direction of the backrest, looking down at the floor. Then do the other arm. Relax. Take some deep breaths. Repeat.
- C. Clasp your hands together on your stomach. Turn your head as far as possible towards each side, stretching the neck muscles. Relax. Take some deep breaths. Repeat.

The most important thing you do on the table is to relax, to let your muscles and tendons stretch. The muscles let “go” of their tension little by little as you “hang” and relax, and by this you regain your flexibility. There is also a beneficial “negative” pressure on your discs, a “suction” action, achieved from just “hanging”/relaxing in 15 degree.

We recommend that you use Mastercare Back-A-Traction one to three times per day, two to fifteen minutes per session.

It is important that you relax and take it easy for a few minutes after each session. You may find it very beneficial just to invert and relax during an entire session, and also to use Back-A-Traction face down/in prone position. In many low back disc problem situations, it is very beneficial to use the table face down/in prone position, especially just before going to bed. We also want you to be prepared that you may experience “popping” in your back as your spine is being aligned.

OBS! Prone position/face down is ONLY pertaining for Models A – Not B or Mini.

Feel free to contact us if you have any questions.

Please note:

For commercial use (hands-on treatments) of Model CN-A1, CN-A2, CN-A3, CN-AH6, it is necessary to utilize the accessory Support Leg, SL1 for the warranty to be valid.

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Exclusive Limited Warranty

Swedish Backcare System, Inc., warrants to the purchaser that should this product be defective by reason of improper workmanship and/or material, warranty repair will be provided without charge for 90 (ninety) days from the date of the original retail purchase, and we will repair or replace, at our option, any defective parts without charge for 2 (two) years from the date of the original retail purchase. Purchaser must notify Swedish Backcare System, Inc., of any defect within the applicable warranty time period in order to qualify for warranty service.

ANY EXPRESS WARRANTY NOT PROVIDED HEREIN, AND ANY REMEDY WHICH, BUT FOR THIS PROVISION, MIGHT ARISE BY IMPLICATION OR OPERATION OF LAW, IS HEREBY EXCLUDED AND DISCLAIMED. UNDER NO CIRCUMSTANCE SHALL SWEDISH BACKCARE SYSTEM, INC., BE LIABLE TO PURCHASER OR ANY OTHER PERSON FOR ANY INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES, WHETHER ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, OR OTHERWISE.

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WARRANTY REGISTRATION CARD

The warranty registration card should have the information necessary to identify the original purchaser, the product type, and identification number. Accordingly, the following information should be included:

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Please return to:

**Swedish Backcare System, Inc.
4055 Pointe Vecchio Circle
Cumming, GA 30040**

**Phone: 770-888-9796
Fax: 770-888-9797**

Model and Serial No. found on the Square-tube beneath the backrest